

# COPY DENTURES

new dentures made by reproducing the features of an existing set

---

## Applications:

- Satisfaction with existing denture
  - Need for fewer clinical visits and treatment time
  - Need for minor modifications
  - Need for a spare denture
  - Material degradation
  - Loss of retention
- 

## Workflow:

- Functional wash impression in occlusion with coeflex (if needed)
  - Trim impression material on buccal side of teeth
- 3D scan 360 degrees
- Send to lab for design/production
- Final placement in patient's mouth

**Turnaround: 2-3 days in lab**

# COPY DENTURES

new dentures made by reproducing the features of an existing set

---

## Applications:

- Satisfaction with existing denture
  - Need for fewer clinical visits and treatment time
  - Need for minor modifications
  - Need for a spare denture
  - Material degradation
  - Loss of retention
- 

## Workflow:

- Functional wash impression in occlusion with coeflex (if needed)
  - Trim impression material on buccal side of teeth
- 3D scan 360 degrees
- Send to lab for design/production
- Final placement in patient's mouth

**Turnaround: 2-3 days in lab**